



EVac[®] 70 Xtra Plasma Wands

Tonsillectomy and Adenoidectomy with One Device

Versatile surgical instruments for the removal of both tonsils and adenoids.

EVac 70 Xtra Plasma wands feature:

Triple-active electrode configuration

- Removes tissue for both tonsillectomy and adenoidectomy


Integrated saline and suction port

- Allows for quick and easy OR set-up

Longer and more malleable shaft

- Improved access to the choanae during adenoidectomy

Enhanced electrode material extends Wand life by approximately 30%¹ (only on the EVac 70 Xtra HP)



System	Catalog No.
Coblator [®] II	EIC8000-01
EVac 70 Xtra Wand	EIC5872-01
EVac 70 Xtra HP Wand	EIC5874-01
EVac Bending Tool	E2020-01

¹ Test Report P/N 18391

ArthroCare ENT
 680 Vaqueros Avenue
 Sunnyvale, CA 94085-3523
 main phone 408-736-0224
 order entry **800-797-6520**
 order entry fax 888-994-2782
www.arthrocareENT.com

ArthroCare Europe AB
 Skeppsbron 2
 111 30 Stockholm
 Sweden
 phone +46 8 546 172 00
 fax +46 8 546 172 39
info@arthrocare.se

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

© 2010 ArthroCare Corporation. All rights reserved. All trademarks, service marks, and logos listed herein are registered and/or unregistered trademarks of ArthroCare Corporation, its affiliates, subsidiaries, or a third party who has licensed its trademarks to ArthroCare. ArthroCare ENT, ArthroCare Sports Medicine, and ArthroCare Spine are affiliates of ArthroCare Medical Corporation.