



PROcise® Plasma Wands

Versatility for Tonsillectomy and Adenoidectomy

Tonsillectomy and adenoidectomy with one device.

Advanced Features Include:

Robust active electrode configuration

- Removes tissue for both tonsillectomy and adenoidectomy
- The PROcise XP features triple-wire active electrodes
- The PROcise EZ features a flat screen active electrode

active triple-wire electrodes



PROcise XP

Decreased diameter on working shaft of Wand

- Increased surgical field visualization—the PROcise EZ features the smallest shaft diameter among all ENT Coblator® wands
- Same working length as EVac® wands
- Bendable shafts for easy access to choanae during adenoidectomy

active flat screen electrode



PROcise EZ

Unique saline delivery system

- Provides tip with saline regardless of orientation
- Unobstructed port clears surgical field to reduce clogging



System	Catalog No.
Coblator® II	EC8000-01
PROcise EZ	EIC8870-01
PROcise XP	EIC8872-01

ArthroCare ENT
 680 Vaqueros Avenue
 Sunnyvale, CA 94085-3523
 main phone 408-736-0224
 order entry **800-797-6520**
 order entry fax 888-994-2782
www.arthrocareENT.com

ArthroCare Europe AB
 Skeppsbron 2
 111 30 Stockholm
 Sweden
 phone +46 8 546 172 00
 fax +46 8 546 172 39
info@arthrocare.se

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

© 2010 ArthroCare Corporation. All rights reserved. All trademarks, service marks, and logos listed herein are registered and/or unregistered trademarks of ArthroCare Corporation, its affiliates, subsidiaries, or a third party who has licensed its trademarks to ArthroCare. ArthroCare ENT, ArthroCare Sports Medicine, and ArthroCare Spine are affiliates of ArthroCare Medical Corporation.