



# PROcise™ EZ View Sinus Wand with integrated ablation, suction, and bipolar hemostasis

*Bringing Enhanced Precision and Visibility to Sinus Surgery*

A single, versatile instrument for your sinus surgery procedures, the PROcise EZ View Sinus Wand brings you precise soft tissue removal offering optimal surgeon visibility and controlled hemostasis in a dry field.

## Features of the PROcise EZ View Sinus Wand include:

### Triple Active Electrode Configuration

- Integrated ablation and suction for controlled soft tissue removal
- Connects to most navigation tracking systems



### Integrated Bipolar Low Temperature Hemostasis

- Variable bipolar coagulation for optimal hemostatic control

### Decreased Diameter on Working Shaft of Wand

- Increased shaft malleability
  - Provides surgeon flexibility during procedures
- Optimal surgeon visibility provides safe navigation and orientation
- Clear access to target tissue in a virtually bloodless operating field
  - Helps maintain orientation in complicated FESS cases for improved patient safety

### Unique Saline Delivery System

- Gravity defying technology ensures the tip is receiving saline flow regardless of orientation
- Strong suction to keep surgical field clear



| Product             | Catalog No. |
|---------------------|-------------|
| Coblator® II System | EC 8000-01  |
| PROcise EZ View     | EIC 8875-01 |

**ArthroCare ENT**  
 680 Vaqueros Avenue  
 Sunnyvale, CA 94085-3523  
 order entry **800-797-6520**  
 phone 408-736-0224  
 order entry fax 888-994-2782  
[www.arthrocareENT.com](http://www.arthrocareENT.com)

**ArthroCare Europe AB**  
 Baggensgatan 25  
 111 31 Stockholm  
 Sweden  
 phone +46 8 546 172 00  
 fax +46 8 546 172 39  
[info@arthrocare.se](mailto:info@arthrocare.se)

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

© 2008 ArthroCare Corporation. ArthroCare, PROcise, Coblator and Coblation are trademarks and/or registered trademarks of ArthroCare Corporation.