



# EVac<sup>®</sup> Xtra HP Plasma Wand with Integrated Cable



*Enhanced Durability and Coagulation*



An “extra-strength” version of the popular EVac 70 Xtra Plasma Wand that delivers extended electrode life and sustained coagulation:

### Triple active electrode configuration

- Tungsten triple-wire electrode extends useful Wand life by 30%
- Maintains thermal penetration in COAG mode throughout Wand life
- Integrated saline line for quick and easy set-up
- Malleable shaft for maximum access



Product	Catalog No.
Coblator II System	EC 8000-01
EVac Xtra HP	EIC 5874-01
EVac Bending Tool	E 2020-01

**ArthroCare ENT**  
680 Vaqueros Avenue  
Sunnyvale, CA 94085-3523  
order entry **800-797-6520**  
phone 408-736-0224  
order entry fax 888-994-2782  
[www.arthrocareENT.com](http://www.arthrocareENT.com)

**ArthroCare Europe AB**  
Baggensgatan 25  
111 31 Stockholm  
Sweden  
phone +46 8 546 172 00  
fax +46 8 546 172 39  
info@arthrocare.se

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.  
© 2008 ArthroCare Corporation. ArthroCare, EVac, Coblator and Coblation are trademarks and/or registered trademarks of ArthroCare Corporation.