



## EVac® 70 Xtra Plasma Wand with Integrated Cable

*One versatile device for Tonsillectomy  
and Adenoidectomy procedures*

A single, versatile surgical instrument for the effective removal of both tonsils and adenoids that provides your patients with minimal pain and fast recovery.

### EVac 70 Xtra Plasma Wand features include:

#### Triple active electrode configuration

- Easily performs tonsillectomy and adenoidectomy
- Variable coagulation for optimal hemostasis
- Integrated suction port

#### Longer and more malleable shaft

- Improved access to the choanae during adenoidectomy
- Increased malleability and compatibility with the EVac Bending Tool

#### Integrated saline line

- Allows for quick and easy OR set-up
- 10' long in length
- Delivers the correct amount of saline to the Wand
- No need to search for additional IV tubine—plug the integrated IV line directly into the saline bag for immediate saline delivery



System	Catalog No.
Coblator® II	EC8000-01
EVac 70 Xtra Wand	EIC5872-01
EVac Bending Tool	E2020-01

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.  
© 2007 ArthroCare Corporation. ArthroCare, EVac, Coblator and Coblation are trademarks and/or registered trademarks of ArthroCare Corporation.

**ArthroCare ENT**  
680 Vaqueros Avenue  
Sunnyvale, CA 94085-3523  
order entry **800-797-6520**  
phone 408-736-0224  
order entry fax 888-994-2782  
[www.arthrocareENT.com](http://www.arthrocareENT.com)

**ArthroCare Europe AB**  
Baggensgatan 25  
111 31 Stockholm  
Sweden  
phone +46 8 546 172 00  
fax +46 8 546 172 39  
[info@arthrocare.se](mailto:info@arthrocare.se)