



## PROcise™ EZ Plasma Wand with Integrated Cable and Suction Line

*Bringing Enhanced Precision and Access to Coblation®*



A single, versatile instrument for your tonsillectomy and adenoidectomy procedures, the PROcise EZ Wand brings you the next generation of ArthroCare ENT Plasma Wands.

### Advanced Features of the PROcise EZ Wand include:

#### Robust active electrode configuration

- Effectively removes tissue for both tonsillectomy and adenoidectomy
- Variable coagulation for optimal hemostasis


#### Decreased diameter on working shaft of Wand

- Increased surgical field visualization – ideal for younger patients
- Same working length as EVac® Wands
- Increased shaft malleability – allows for greater surgeon flexibility during procedures

#### Unique saline delivery system

- Gravity defying technology ensures the tip is receiving saline flow regardless of orientation
- Larger, stronger suction port than the EVac Wands
- Unobstructed port efficiently clears surgical field to reduce clogging

#### Integrated 10' saline line with spike and integrated cable for OR convenience and set-up



System	Catalog No.
Coblator® II	EC8000-01
PROcise EZ	EIC8870-01



CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.  
 © 2006 ArthroCare Corporation. ArthroCare, Coblator, Coblation, EVac, and PROcise EZ are trademarks and/or registered trademarks of ArthroCare Corporation.

**ArthroCare ENT**  
 680 Vaqueros Avenue  
 Sunnyvale, CA 94085-3523  
 order entry **800-797-6520**  
 phone 408-736-0224  
 order entry fax 888-994-2782  
[www.arthrocareENT.com](http://www.arthrocareENT.com)

**ArthroCare Europe AB**  
 Baggensgatan 25  
 111 31 Stockholm  
 Sweden  
 phone +46 8 546 172 00  
 fax +46 8 546 172 39  
[info@arthrocare.se](mailto:info@arthrocare.se)