




PROcise™ Plasma Wands with Integrated Cable and Suction Line

Bringing Enhanced Precision and Access to Coblation®



A single, versatile instrument for your tonsillectomy and adenoidectomy procedures, the PROcise Wands bring you the next generation of ArthroCare ENT Plasma Wands.

Advanced Features of the PROcise Wands include:

Robust active electrode configuration

- Effectively removes tissue for both tonsillectomy and adenoidectomy
- Variable coagulation for optimal hemostasis



PROcise EZ

Decreased diameter on working shaft of Wand

- Increased surgical field visualization – ideal for younger patients
- Same working length as EVac® Wands
- Increased shaft malleability – allows for greater surgeon flexibility during procedures



PROcise XP

Unique saline delivery system

- Gravity defying technology ensures the tip is receiving saline flow regardless of orientation
- Larger, stronger suction port than the EVac Wands
- Unobstructed port efficiently clears surgical field to reduce clogging

Integrated 10' saline line with spike and integrated cable for OR convenience and set-up

System	Catalog No.
Coblator® II	EC8000-01
PROcise EZ	EIC8870-01
PROcise XP	EIC8872-01

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.
 © 2007 ArthroCare Corporation. ArthroCare, Coblator, Coblation, EVac, PROcise EZ and PROcise XP are trademarks and/or registered trademarks of ArthroCare Corporation.

ArthroCare ENT
 680 Vaqueros Avenue
 Sunnyvale, CA 94085-3523
 order entry **800-797-6520**
 phone 408-736-0224
 order entry fax 888-994-2782
www.arthrocareENT.com

ArthroCare Europe AB
 Baggensgatan 25
 111 31 Stockholm
 Sweden
 phone +46 8 546 172 00
 fax +46 8 546 172 39
info@arthrocare.se

